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**CONTRA-ANGLE SURGICAL INSTRUMENTS**

**PRODUCT IDENTIFICATION**

Surgical instruments related to implant systems manufactured by Kristal Srl are medical devices intended for use in the oral cavity, for temporary use (continuous duration not exceeding 60 minutes), reusable, in NOT STERILE packaging.

The functions of the surgical instruments are:

- preparation of sites for Bio Implant implants
  - the insertion of the plants in the sites
  - the tightening and unscrewing of all connection screws (cap screws, healing transgingival screws, abutment screws, abutments, prosthetic screws, transfer screws, etc...).
- The surgical instruments manufactured by Kristal Srl are intended for use with dental implants, also manufactured by Kristal Srl. The use of surgical instruments for operations with implants other than those manufactured by Kristal Srl limits the responsibility of Kristal Srl and cancel the product warranty (see section "Defective Product Liability and Warranty Terms" below). We are not responsible for the use of non-original equipment.

**DESCRIPTION AND USE**

The information in these operating instructions complete the instructions in the catalogs/manuals. Request a copy from Kristal Srl if you do not have it. All the devices are identified by the instrument code, which is reported with a laser marking on the body of the same. If the space does not allow to return the complete code, the elements that allow to identify the device in a univocal way (eg diameter or length) are however marked.

Within the implant systems Bio Implant is defined a color code that helps to identify the diameters of the instruments to be used depending on the implant diameter or the measurement of the platform. The color code legend is explained in the catalogs and surgical manuals of the individual implant systems. The instruments covered by these operating instructions (surgical drills and extension cutters, drivers, surgical trays) are intended for mechanical use, ie they have a shank with a contra-angle attachment and must be used with a suitable micromotor. Incorrect insertion of the instruments into the handpieces can lead to instrument vibrations, excessive rotations, premature wear and shank bending. It is recommended to use only surgical micromotors suitable for use. It is recommended to have the micromotors periodically checked by the manufacturers, according to the individual indications of the same, to prevent possible malfunctions (eg axis displacements of the transmission shafts, worn or malfunctioning pliers, etc...). Failure to follow the instructions provided may cause surgical problems and damage to the patient's health.

**DRILLS AND EXTENSION CUTTERS**

They are sharp instruments, to be used for the preparation of surgical sites for dental implants. They have different morphologies: lanceolate, super-cut, conical, bone trap, countersinkers, bone profilers. Depending on the implant system to which they belong, they are or are not equipped with internal irrigation, have notches of depth to allow the determination of the working depth, have or not the predisposition to hook a stop depth. Refer to the catalogs of the individual implant systems for detailed technical specifications. The drills produce a hole that is longer than the insertion depth of the instrument. This greater length is linked to the measurement of the tip of the instrument. The exact measurements of the depth of the holes determined by the drills are shown in all the catalogs and surgical manuals of the implant systems.

It is recommended to use the rotation speeds indicated in the single catalogs or surgical manuals, to avoid the development of bone necrosis. Lever movements increase the risk of fracture of instruments, therefore they must be avoided. In general, sudden changes in speed must be avoided. Never apply pressure such as to stop the rotation of the instrument by force. This could lead to excessive heat increase in the tissues affected by the cut, resulting in bone necrosis, and ruining both the instrument and the device used (micromotor). This could also mean breaking the instrument itself. It is also recommended to work intermittently, to avoid overheating and wear of the working

part and undue increase of heat in the tissues affected by the cut. The use of appropriate refrigeration liquid is recommended. In the absence of adequate irrigation bone necrosis may be generated.

The consumption of the drills depends largely on the type and density of the milled bone: harder bone means greater wear of the instruments. It is recommended, for greater safety and prudence with respect to the wear resistance of the device, that the drills are used for no more than 20 cycles of work or earlier if the tools lose their cutting capacity. The recommended 20 cycles represent an average figure. It is recommended to check the maintenance status of the residual cutting capacity after each intervention.

Kristal Srl assumes no responsibility in case of excess use. The drills must never be resharpened before use. Never use damaged, bent or worn tools. An extension cutters are available, to be used when the total length of the instruments is too short depending on the presence of adjacent teeth that do not allow the head to pass through the handpiece. If you are using it, make sure that the shank of the cutter is inserted correctly and completely. An incomplete insertion determines eccentricity of rotation of the cutter.

**STOP FOR DRILLS**

They are devices to be inserted on the drills prepared to receive them. They allow to limit the working length of a cutter at a predetermined height. The measurements are shown on the catalogs, on the surgical manuals and screen printed in the surgical trays. Always check that the stop is inserted at the desired height. An incomplete insertion can reduce the height of the preparation. Possible insertion difficulties can be solved by slightly loosening the stop tabs, using tweezers. Also verify that the stop performs sufficient retention. Low weak retention can cause the instrument to fall off the drill during surgery. In the event that the stops diminish the retention capacity, simply tighten the fins slightly, either manually or using tweezers.

**DRIVER**

They are devices with the dual function of carrier, which allow you to take the implants from their packaging without contaminating them, that is without touching the surface, and to transport them to the oral cavity without ever touching them, and the function of screwing driver, transmitting the rotary motion from the micromotor to the implants allowing the screwing in the prepared sites.

Lever movements should be avoided because they increase the risk of fracture. The drivers are different, depending on the different implant systems. The technical details of each individual system can be found in the surgical manuals and catalogs. It is recommended to read it before use.

Screwdrivers that are used for screwing abutment screws or prosthetic screws must be used with tightening torque control:

- Through screws to tighten on the pillar and abutment implants: 20-25 Ncm
- Through screws for tightening prosthetic superstructures to abutments: 20-25 Ncm

- Screwing of components that screw directly to the implants (eg spherical attachments, some types of abutments that do not have a through screw but make a single body with the screw): 30 Ncm

- Through screws for clamping superstructures directly on the implants (without using intermediate abutments): 20-25 Ncm

In the case of screwing of transmucous healing screws it is recommended not to exceed 8-10 Nm of torque. Too much tightening torque can weaken the mechanical structure of the screws and compromise prosthetic stability, with possible damage to the implant connection.

**SURGICAL TRAY**

They are practical trays within which all the necessary instruments for the surgery and the prosthesis of Bio Implant dental implants are available and organized in an ergonomic way.

Please refer to the individual catalogs and user manuals for detailed information about the contents of the individual kits and the procedures for use. For the purposes of Directive 93/42/EEC implemented in Italy with D.L. 46/97 of 26/3/97, Annex I and related updates and amendments, Kristal Srl declares the manufacturer of the surgical instruments for Bio Implant dental implants and identifies the risk class as follows:

- Surgical drills (lanceolate, super-cut, conical, bone pickers, countersinkers, bone profiling) and Drill extension, Stop for drills, Drivers, and Screwdrivers/Screwdrivers: Surgical medical devices invasive for temporary use (for a duration of less than 60 continuous minutes), reusable, in NOT STERILE packaging, Risk Class 2a according to the rules of Annex IX
- Surgical Tray: Medical Devices, in NOT STERILE packaging, reusable; Hazard Class 2 because they contain part or all of the aforementioned instruments, necessary for the surgery and the prosthesis of Bio Implant dental implants. The use and handling of the product is reserved for medical and dental personnel with the necessary qualifications and professional preparation.

**IDENTIFICATION OF THE MANUFACTURER**

The manufacturer of surgical instruments for dental implants covered by these Instructions for Use is: Kristal Srl  
Registered office: Via Cogorano, 25 - 57123 Livorno (LI) - Italy  
Operating office: Via Edison, 14 - 20090 Trezzano S/N (MI) - Italy  
Tel. +39 024453401 - info@kristalrsl.net - www.kristalrsl.net

**RAW MATERIAL USED**

The materials used for the production of surgical instruments for dental implants manufactured by Kristal Srl have been selected on the basis of the properties indicated for their intended use, in compliance with Directive 93/42/EEC implemented in Italy by Law 46/97, Annex I Essential Requirements, point 7.1. They are produced, depending on the type of component, in:

- Titanium grade 5 (Ti6Al-4V)
- AISI 630 steel

It is recommended to check with patients any allergy to the listed raw materials. Please refer to the website www.kristalrsl.net for detailed technical data sheets of all the materials used, for the verification of the relative chemical compositions and for the physical-mechanical characteristics.

**WARNINGS**

The Bio Implant surgical instruments are sold in NOT STERILE packaging. Be-fre using them, they must be cleaned, disinfected and sterilized according to the instructions given below. Failure to comply with this warning may cause patient infections. It is advisable to collect and file a complete clinical, radio-logical and radiographic documentation.

The code, the description of the content and the lot number are shown on each package. These same data are also indicated on the labels placed inside the packs and will always be quoted by the doctor for any communication on the matter.

When handling devices, both during use and during cleaning and sterilization operations, it is recommended to always use surgical gloves for personal protection against bacterial contamination. Failure to follow these rules can result in cross infections. The package complies with European standards.

**CONTRAINDICATIONS**

In evaluating the patient, in addition to considering the suitability for implant prosthetic rehabilitation, it is generally necessary to take into account the contraindications valid for dental surgery.

- Among these are mentioned:
- Alterations of the blood coagulation chain, therapies performed with anticoagulants
  - Scarring or bone regeneration disorders such as:  
- Non-compensated diabetes mellitus
  - Metabolic or systemic diseases of the part that compromise the regeneration of tissues with particular impact on cicatrization and bone regeneration
  - Abuse of alcohol and tobacco and use of drugs
  - Immunosuppressive therapies such as chemotherapy and radiotherapy
  - Infections and inflammations such as: periodontitis, gingivitis
  - Lack of oral hygiene
  - Inadequate motivation
  - Defects of occlusion and/or joint as well as insufficient interocclusal space
  - Inadequate abutment/prosthetic
- The drills must be used with caution in cases of low bone density, and the

sites must be adequately underprepared. Preferably use bone condensers. The insertion of implants and implant prostheses is contraindicated in patients with poor general health, poor or insufficient oral hygiene, impossibility or lack of control over general conditions, or who have previously undergone organ transplants.

In addition, psycholabile patients, or those who abuse alcohol or drugs, with insufficient motivation or insufficient cooperation must be discarded. Patients with bad periodontal status must be previously treated and recovered. In case of lack of bone substance or poor quality of the receiving bone, such that the stability of the implant can be impaired, an appropriate guided tissue regeneration must be carried out beforehand. They also represent a contraindication: titanium allergy, acute or chronic infectious diseases, chronic subacute maxillary osteitis, systemic diseases, endocrine disorders, diseases with consequent microvascular disorders, pregnancy, lactation, previous exposure to radiation, hemophilia, granulocytopenia, steroid use, diabetes mellitus, renal failure, fibrous dysplasia. Implants intended to support the prosthesis are medical devices that are inserted into the oral cavity as part of a surgical procedure, and as such provide for further restrictions of use, for which reference should be made to the Instructions for Use relating to implant fixtures.

**COLLATERAL SYMPTOMS**

Among the manifestations that accompany the surgery can occur:

- Temporary local swelling, edema, hematomas.
- Temporary limitations of sensitivity.
- Temporary limitations of masticatory functions.
- Post-operative microhaemorrhages in the next 12/24 hours.

**CLEANING/DISINFECTION/STERILIZATION/CONSERVATION**

All Bio Implant dental implants are sold in NON-STERILE condition. Before use, they must be cleaned, disinfected and sterilized. These processes must be performed before the first use, and before each subsequent reuse. Repetition of the processes described in this section has a minimal effect on these devices. The correct functionality of the instruments must always be checked before their use. In the presence of signs of wear, the instruments must be replaced immediately with new devices. In particular, it is recommended to always check the correct retention of the screwdrivers inside the mounting wells on the head of the screws that must be removed and screwed with them. Failure to follow these indications may result in cross infection and intraoperative complications.

**CLEANING**

Containers and transport to be used for washing: there are no particular requirements.

In the case of automated cleaning: use an ultrasonic tank using a suitable cleaning solution. It is recommended to use only neutral detergents. The concentration of the solution and the duration of the washing must follow the manufacturer's instructions. Use demineralized water to prevent the formation of stains and stains.

When unloading, check the recesses of the devices, the holes, etc. to verify the complete removal of residues. If necessary, repeat the cycle or use manual cleaning.

In the case of manual cleaning: use a suitable neutral detergent, following the manufacturer's instructions for use. Brush the products with soft brushes, under plenty of running water. Using the brush, apply the cleaning solution to all surfaces. Rinse with distilled water for at least 4 minutes. Make sure that running water passes in abundance through any holes.

In the case of drills with internal irrigation, use the special pins supplied with the handpieces to make sure that the irrigation holes have been completely cleaned and free from the residues of bone fistules or biological tissues. After rinsing, completely dry the devices and bag them in suitable sterilization bags. If a drying cycle is carried out as part of the cycle of a washing and disinfecting appliance, do not exceed 120° C.

**STERILIZATION**

In a vacuum autoclave, and sterilize in the following way:

- Temperature = 121 ÷ 124° C, with a minimum autoclave cycle of 20 minutes and a drying cycle of 15 minutes.

**STORAGE**

After sterilization, the product must remain in the bags used for sterilization. Envelopes should only be opened immediately before re-use. The sterilization pouches are normally able to maintain sterility inside them, except for damage to the envelope. Therefore, be careful not to use components if the envelopes in which they were stored are damaged and re-sterilize them in new bags before re-use. The period of storage of the sterilized products inside the envelopes must not exceed that recommended by the manufacturer of the envelopes themselves. The product must be stored in a cool and dry place, away from direct sunlight, water and sources of heat.

**NONMATIVE REQUIREMENTS**

The design and production of the surgical components is performed in compliance with the most up-to-date harmonized standards and directives regarding the materials used, the production processes, the information supplied and the packaging.

**DISPOSAL PROCEDURES**

Surgical instruments, if used, must be assimilated for their disposal to biological waste, according to local regulations.

**RESPONSIBILITY OF THE DEFECTIVE PRODUCT AND GUARANTEE TERMS**

Optimal patient care and attention to their needs are necessary conditions for implant success and it is therefore necessary to carefully select the patient, inform him of the inherent risks and duties associated with the treatment and encourage him/her to cooperate with the dentist for the successful outcome of the treatment itself.

Therefore, the patient must maintain good hygiene, confirmed during check-ups and check-ups; it must always be ensured and documented as, moreover, the indications and instructions of the doctor must be observed and documented.

The instructions provided by Kristal Srl. they are available at the time of treatment and accepted by the dental practice; it is necessary to observe and apply them at all stages of treatment: from the patient's medical history to post-operative check-ups.

The warranty covers only the established production defects, after sending the piece identified by item code and lot, within the warranty period. The warranty terms are available on the website www.kristalsrl.net

**DATE AND VALIDITY OF THESE INSTRUCTIONS FOR USE**

These instructions for use are valid and effective from January 2019.

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