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## INSTRUCTIONS FOR USE OF BIO IMPLANT PROSTHETIC COMPONENTS

### Description

The prosthetic components are made of Titanium C.P. grade 4 or grade 5 (Ti6Al4V) that can be subjected to anodizing, gold alloy, acrylic castable, biocompatible polymer (PEEK) and ceramic (Zirconia).

The brand Bio Implant prosthetic components are dedicated to the reference implant line (Plus, K-Plus, Core, K-Core, Core V2, K-core V2). For the V2 line characterized by unified prosthetic platform are three different emergency profiles available: Narrow, Regular, Wide; to the exclusion of  $\varnothing 2,9$ , for which is scheduled dedicated prosthetic components.

They give some technical characteristics of the prosthetic made of gold alloy (Platinor N) and biocompatible polymeric material (PEEK).

- prosthetic articles gold alloy
- Chemical composition % (m / m) with a certain emission spectrometer (Spectro Lab S, Spectro, Kleve): Au 60.0, 24.9 Pt, Pd 15.0, 0.1 Ir.
- Reference standard: UNI EN ISO 22674 Dentistry: metallic materials for fixed and removable restorations and appliances.
- Melting range: 1350°C - 1460°C.
- prosthetic articles polymer material PEEK:
- Temporary employment (max. 30 days);
- Thermal properties: glass transition temperature: 150°C;
- Melting temperature: 341°C;
- Short-term operating temperature (30 days): 300°C.
- prosthetic articles in zirconia ceramic material:
- Hardness of Zirconia: 1200/1400 HV;

- Workability: Workable with diamond tools.

All devices are packaged in double container, the first pack is a hermetically sealed envelope or a blister PET thermally sealed with medical paper in the case of V2 line, the second is a cardboard container in which the immediate packaging is placed. Both packs are labeled with the identification code, batch number and other useful information for the proper use of the product. The product is supplied non-sterile clean.

All prosthetic components brand Bio Implant fulfills the essential requirements of Directive 93/42/EEC and subsequent amendments concerning medical devices.

### Indications

The prosthetic components Bio Implant is indicated in all cases of prosthetic-functional and aesthetic rehabilitation both maxillary and mandibular, following the use of endosseous dental implants of the same brand. The prosthetic components of the Plus and K-Plus range have internal hex; prosthetic components of the Core line and K-Core feature external hex as well as the V2 line characterized by unified prosthetic platform, with the exception of dental implants  $\varnothing 2,9$  which are provided with prosthetic components dedicated. In particular after plant application we proceed to the closing of the flaps. For osseointegration occurred is carried out before application for the reopening of the healing abutment for the conditioning of the gum and then the only final prosthesis. The prosthesis is made with single crowns, dentures made solid, full or partial dentures attached to plants by appropriate prostheses.

Before use, we recommend sterilization of prosthetic components by performing a standard autoclave sterilization cycle (121°C 1.1 bar 15 min. / 134°C 2.1 bar 4 min.).

### Applications

The prosthetic components must be selected and used depending on the patient's condition, surgical experience and the operative technique. The indications on the use of prosthetic components are listed on the Bio Implant catalogs.

### Contraindications

The contraindications of the prosthetic components of the line Bio Implant Core V2 are linked to those of dental implants of the mark, since it is used later to the system of such devices.

### Warnings

The prosthetic components of the line Bio Implant can only be used by dentists, dental surgeons, maxillofacial surgeons and dental technicians to create dentures. These components are an integral part of a complete system, so it is recommended to use original components.

They should be stored in the original packaging until ready to use in a dry place at room temperature. The device is strictly disposable. Reuse of single use products creates a potential risk of infection in the patient or in those who use it. E' to be affixed on the patient's medical record the identification data of prosthetic parts used. It is advisable to also make, radiological documentation, clinical and complete statistics.

The patient must be informed of the importance of cleanliness and maintenance of the prosthesis. It must then be instructed on how to proper hygiene and maintenance of the prosthesis.

By the holders of oral implants in Titanium they must not be used toothpastes or mouthwashes that contain free fluoride.

In cases of necessity of therapeutic radiation to the head and neck, metal dentures must be removed from the mouth. Accessories for endosseous implantology must be regarded as a special waste and this is taken into account in the event of disposal. And 'the responsibility of ensuring the correct method of disposal.

To obtain the replacement of the goods when the quality does not match the standards of the manufacturer, if the wrong order or if visibly damaged packaging, it is necessary to warn by telephone or by e-mail within 10 days of receipt of goods; making the return of the product in its original filling a regular transport document of the commodity and the authorization request document that comes from your target agent.

### Known characteristics and technical factors known to the manufacturer that could cause danger

Incorrect connection between implant and abutment, prosthetic load will not physiological, incorrect mounting of the cap on the pier with a consequent ingestion of the same, absence of healing cap when necessary.

### Precautions

Must comply with the post-operative procedures and the healing time recommended in the Surgical Manual prior to replacement.

### Side effects

The secondary effects of the components of the prosthetic line Bio Implant are linked to those of the dental implants of the same mark, since it is used following implantation of such devices.

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