# PRIMARY HEALING IMPLANT **PHI**

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## KRISTAL

IMPLANTOLOGY - ORTHODONTICS

Kristal S.r.l.  
Sede legale: via Cogorano, 25 - 57123 Livorno (LI) Italia  
Sede operativa: via Edison, 14 - 20090 Trezzano S/N (MI) Italia  
Tel. +39 024453401 - info@kristal.srl.net - www.kristal.srl.net

- 1 ISTRUZIONI PER L'USO DEGLI IMPIANTI LINEA PHI
  - 2 INSTRUCTIONS FOR USE OF PHI DENTAL IMPLANTS
  - 3 MODE D'EMPLOI DES IMPLANTS DENTAIRE DE LA LIGNE PHI
  - 4 BEDIENTUNGSANLEITUNG FÜR DIE ZAHNIMPLANTATE LINE PHI
  - 5 INSTRUCCIONES PARA EL USO DE LOS IMPLANTES PHI
- 1 ISTRUZIONI PER L'USO DEGLI IMPIANTI LINEA PHI

### Descrizione

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### 2 INSTRUCTIONS FOR USE OF PHI DENTAL IMPLANTS

### Description

PHI CLASSIC and PHI EVO dental implants are made of Grade 5 Titanium according to the ASTM F136. They have an internal hex connection, headquarters for prosthetic components dedicated to specific line.

Each implant is packaged in a double aseptic vial. The primary container is a transparent PE bottle with screw cap with medical silicone insert; inside a titanium cylinder supports the system that exposes only by the neck and its hexagon. The implant collection is carried out using manual key or mount contra handpiece. The cover screw is housed in insert in the cap silicon. All the bottles are placed in a secondary package or a blister PET thermally sealed with medical paper, on which all the main features of the device is affixed bearing a label such as: identification code, description, lot, expiration date and the related symbols.

The dental implants are supplied in a sterile package. Sterilization is performed by Gamma irradiation (25 kGy). Dental implants to mark PHI meet the essential requirements of Directive 93/42/EEC following changes/additions, relatives on medical devices.

### Use of the implants

The PHI dental implants are indicated in all cases of prosthetic-functional and aesthetic rehabilitation, both maxillary and mandibular. The use of the implant system is provided for treatment total or partial edentulism. The CLASSIC line is used for the transmucosal technique; The EVO line is used for the submerged technique.

The duration of the installations also depends on proper maintenance by the patient, which must be fully informed about the

tasks needed. Also it is the obligation of the physician to make periodic written reminders for inspection and maintenance of implants.

The CLASSIC and EVO implants should be considered an integral part a complete implant system that includes PHI in addition to accessories and the instruments also the prosthetic components available in various configurations characterized in the case of the line EVO from three different emergence profiles: high, low or totally closed (no platform-switching).

### Table for minimum implant site

Position	CLASSIC	EVO
17-27	4,5/10	4,5/8
16-26	4,5/13	4,5/10
15-25	4,5/10	4,5/8
14-24	4,0/10	4,0/8
13-23	4,0/13	4,0/10
12-22	3,5/10	3,5/8
11-21	4,0/10	4,0/8
47-37	4,5/10	4,5/8
46-36	4,5/13	4,5/10
45-35	4,5/10	4,5/8
44-34	4,0/10	4,0/8
43-33	4,0/13	4,0/10
42-32	3,5/10	3,5/8
41-31	3,5/10	3,5/8

Caution when using smaller plants are obliged to warn the patient that these systems are out of protocol and their duration may be less than normal (over 80% in ten years).

### Applications

The implants must be selected according to the patient's condition and the surgical technique. Specifically, the surgical technique is explained in the surgical protocol. A high degree of accuracy is required at all times in the preparation of the implant site by respecting the sequence of the various steps in the surgical protocol. The heat trauma caused by friction of the drills must be kept to a minimum by using cutters sharp, in the right consequential way, the appropriate number of round and irrigation with sterile solutions and suitable healing times.

### Contraindications

Insufficient amount of bone tissue or gum or insufficient or inadequate bone quality. Allergy to titanium or its alloys with respect to patient selection to observe the same precautions and contraindications for dental and general surgery. Some of them include: reduced blood clotting (anticoagulants to antiplatelet therapies or, hereditary or acquired coagulation disorders) infection and inflammation of the mouth and the surrounding areas such as periodontitis and gingivitis (acute and chronic infections of the affected area, or infection inflammation of the soft tissues, bone infections, osteomyelitis). More over severe metabolic diseases (decompensated diabetes or inadequately treated, the calcium metabolism disorders, treatment with steroids and other drugs that interfere in bone metabolism, therapy with immune-suppressants, chemotherapy and radiotherapy, active endocrinological diseases, bone diseases, anatomical proximity or interference with nerve routes. Also insufficient gingival soft tissue for recasting, occlusion and / or dental articulation incorrect, limited interocclusal space. Further more psychic disorders, pain syndromes, poor oral hygiene, poor patient motivation to this type of treatment. relative contraindications apply to bruxism, allergic patients, patients with addictions to alcohol, drugs, nicotine.

### Warning

The PHI dental implants may only be used by Dentists Doctors, dentists and maxillofacial surgeons, properly trained in the use of implant technique PHI. Failure to observe the surgical protocol could determine the implant failure. PHI CLASSIC and PHI EVO dental implants are considered an integral part of a complete implant system. Kristal recommend the use of accessories, tools and original prosthetic components. Dental implants must be stored in its original packaging in a dry place, away from heat sources and at room temperature until the time of placement. Dental implants should not be used after the expiration date on the package.

Do not use the product if the package is opened or damaged. The lack of labeling and instructions may result in the use of the product beyond the expiration date and may be a necessary traceability. To modify the specification of the product packaging configuration could compromise sterility and functionality of the device.

Dental implants are sold in sterile and single-use packaging, the reuse of single-use products creates a potential risk of infection in the patient or in those who use it.

The manufacturer is not liable in case of re-sterilization performed by third parties. Reprocessing may compromise the structural integrity of the product or cause it to malfunction. The removal of the plant from its packaging needs of the ideal about the aseptic standards. Contamination of the device may result in serious infection of the patient.

The application of the final prosthesis should occur no earlier than 12 weeks after implant placement in standard conditions.

Kristal recommended to draw up appropriate clinical documentation. To ensure traceability of each plant is necessary to store information such as part number and lot.

Kristal provides a special label that shows the main features of the dental implant inserted to be placed on the patient's medical record.

To obtain the replacement of the goods when the quality does not match the standards of the manufacturer, if the wrong order or if visibly damaged packaging, it is necessary to warn by telephone or by e-mail within 10 days of receipt of goods; making the return of the product in its original filling a regular transport document of the commodity and the authorization request document made that comes from your target agent.

The devices removed from the oral cavity for a biological or mechanical failure must be treated for their disposal in biohazard waste. In general, please refer to the applicable regulations locally.

### Precautions

In order to facilitate the surgical planning is absolutely necessary an accurate study investigation prior to treatment: screening, pa-

noramic radiographs and CT.

### Secondary effects

To date Titanium has no known side effects. Despite the high success rate of Titanium implants, there is a risk of failure which is generally linked to the incorrect surgery techniques and disregards for absolute implantology contraindications.

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