PRIMARY  
HEALING  
IMPLANT

**PHI**

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**KRISTAL**  
IMPLANTOLOGY • ORTHODONTICS

Kristal S.r.l. via Edison, 14  
20090 Trezzano S/N (MI) - Italia  
Tel. +39 024453401 - info@kristalrfl.net

- 1 ISTRUZIONI PER L'USO DELLE COMPONENTI PROTESICHE LINEA PHI
- 2 INSTRUCTIONS FOR USE OF PROSTHETIC COMPONENTS PHI LINE
- 3 MODE D'EMPLOI DES COMPOSANTS PROTHÉTIQUE LINE PHI
- 4 HINWEISE FÜR DIE VERWENDUNG VON COMPONENTS PROSTHETIC LINE PHI
- 5 INSTRUCCIONES DE USO DE LOS COMPONENTES PROTÉSICOS PHI LINE

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## INSTRUCTIONS FOR USE OF PROSTHETIC COMPONENTS PHI LINE

### Description

The prosthetic components are made of grade 5 Titanium according to the ASTM F136, except for the healing screws made of Titanium grade 4 CP, and it can be coated in Titanium Nitride (TiN); Acrylic castable and biocompatible polymer (PEEK).

The prosthetic components PHI are dedicated to the implant lines of CLASSIC reference (cemented) and EVO (screwed).

The prosthetic components meet the essential requirements Directive 93/42/EEC and subsequent amendments concerning medical devices.

They give some technical characteristics of the prosthetic made of biocompatible polymeric material (PEEK):

- Temporary employment (max 30 days);
- Thermal properties: glass transition temperature: 50°C;
- Melting temperature: 341°C;
- Short-term operating temperature (30 days): 300°C.

### Packaging

All items are packaged as follows: prosthetic components are housed in a heat-sealed blister with medical cards contained in a cardboard box. Both packs have proper labeling indicating the reference item number, lot number and other useful information for the proper use of the product. The product is supplied non-sterile clean.

## Indications

The prosthetic components PHI is indicated in all cases of prosthetic-functional and aesthetic rehabilitation both maxillary and mandibular, following the use of endosseous dental implants of the same brand. CLASSIC after plant application is part of a closure screw of corresponding dimensions to the system; while for the line EVO proceed to the closing of the flaps and to osseointegration occurred is carried out before application for the reopening of the healing abutment for the conditioning of the gum, and then the only final prosthesis. The prosthetic components of the PHI line are compatible to all installations of the line EVO or CLASSIC to which they are dedicated. The line CLASSIC involves the use of definitive abutments cemented within the implant cavity through the use of cold cement metal/metal, unless solutions are used for bars or overdentures coupling head. The EVO line involves the use of abutments with through screw with three different emergency profiles: high, low, or total closure (no platform-switching); also for this line are various solutions for the realization of prostheses or overdentures on bar coupling head.

The prosthesis is made with single crowns, dentures made solid, full or partial dentures attached to plants by appropriate prosthesis.

Before use, we recommend sterilization of prosthetic components by performing a standard autoclave sterilization cycle (121°C 1.1 bar 15 min. / 134°C 2.1 bar 4 min.).

## Applications

The prosthetic components must be selected and used in function the patient's condition, surgical experience and the operative technique. The indications on the use of prosthetic components PHI are listed on the respective catalogs.

## Contraindications

The contraindications of the PHI line prosthetic components are linked to those of the dental implants of the reference line, since it is used following implantation of such devices.

## Warnings

The prosthetic components of the PHI line can only be used by dentists, dental surgeons, maxillofacial surgeons and dental technician for creating prosthetic devices. These components are an integral part of a complete system, so it is recommended to use original components.

Should be stored in its original packaging until ready to use in a dry place at room temperature. The device is strictly disposable. Reuse of single use products creates a potential risk of infection in the patient or in those who use it. To be affixed on the patient's medical record the identification data of prosthetic parts used. It is advisable to also make, radiological documentation, clinical and complete statistics.

It's available a prosthetic manual supplied on request in which are specified all the clinical and laboratory procedures for the installation of the prosthetic components. The patient must be informed of the importance of cleanliness and maintenance of the prosthesis. It must then be instructed on how to proper hygiene and maintenance of the prosthesis.

From the carriers of oral implants in Titanium must not be used toothpastes or mouthwashes that contain free fluoride.

In cases of need for therapeutic radiation to the head and neck, metal implants should be removed from the mouth.

For replacement of the goods when the quality does not match the standards of the manufacturer, if the wrong order or if visibly damaged packaging, it is necessary to warn by telephone or by e-mail within 10 working days of receipt of goods; making the return of the product in its original filling a regular transport document of the commodity and the authorization request document made that comes from your target agent.

The devices removed from the oral cavity for a biological or mechanical failure must be treated for their disposal in biohazard waste. In general, please refer to the applicable regulations locally.

## Known features and the technical factors known to the manufacturer that could cause danger

Incorrect connection between fixture and abutment; prosthetic loading not physiological; Incorrect mounting of the cap on the pier with a consequent ingestion thereof.

## Precautions

Must comply with the post-operative procedures and the healing time recommended in the Surgical Manual prior to replacement.

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