Kristal S.r.l.
Sede legale: via Cogorano, 25 - 57123 Livorno (LI) Italia
Sede operativa: via Edison, 14 - 20090 Trezzano S/N (MI) Italia
Tel. +39 024453401 - info@kristalimplant.net - www.kristalimplant.net

- 1 ISTRUZIONI PER L'USO DEGLI IMPIANTI LINEA BIO IMPLANT
- 2 INSTRUCTIONS FOR USE OF DENTAL IMPLANTS BIO IMPLANT LINE
- 3 MODE D'EMPLOI DES INSTALLATIONS DE LIGNE BIO IMPLANT
- 4 HINWEISE FÜR DIE VERWENDUNG VON ANLAGEN LINE BIO IMPLANT
- 5 ISTRUCCIONES PARA EL USO DE LOS IMPLANTES DE LA LINEA BIO IMPLANT

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2 INSTRUCTIONS FOR USE OF DENTAL IMPLANTS BIO IMPLANT LINE

Description

Dental implants of Bio Implant line are made of Titanium C.P. grade 4 according to ASTM F67. Shipped with an external hex connection for the Plus and K-Plus line, inside the Core and K-Core line, with unified internal prosthetic platform for the Core V2 and K-Core V2 line, to Ø2.9 implant exclusion. Bio Implants Plus, K-Plus, Core, K-Core, Core V2 Ø2.9 are available with pre-assembled mounting device to the system and can be used in conjunction with dedicated prosthetic components.

Each implant is packaged in a double aseptic vial:

- For the Plus and K-Plus line, Core and K-Core line, the primary container is a transparent PE bottle with screw cap equipped with a label for encoding color corresponding to the diameter of the implant; inside a cylinder Titanium also capable 4, supports the implant mount that only exposes the control panel. The mounting device (mount) is preassembled on the system directly in the clean room, and allows the practitioner to pick up and place the same in conditions of absolute sterility. The plant with the mount is removed from its housing by means of a screwdriver (manual, ratchet and contra). The mount is screwed into the implant through a passing screw to be discarded after insertion. The mount can be used as a temporary mesostructure as realized in Titanium C.P. grade 4. The cover screw is included in the package, and is located laterally to the system.
- For the Core V2 and K-Core V2 line, the primary container is a transparent PE bottle with screw cap with insert colored medical grade silicone corresponding to the implant diameter; inside a cylinder Titanium also capable 4, it supports the implant that exposes only for the neck and the relative hexagon, except for the implant Core V2 Ø2.9 that presents the pre-assembled mounting device. The cover screw is housed in insert in the cap silicon.

All the bottles are placed in a secondary package or a blister PET thermally sealed with medical paper, on which all the main features of the device is affixed bearing a label such as: identification code, description, nr. Lot, expiration date and the related symbols.

Dental implants are supplied sterile. The sterilization is carried out with irradiation by gamma rays (25 kGy). Dental implants brand Bio Implant meet the essential requirements of Directive 93/42/EEC and subsequent amendments on medical devices.

Indications for Use

Bio Implant dental implants are indicated in all cases of prosthetic-functional and aesthetic rehabilitation, both maxillary and mandibular.

The use of the implant system is provided for treatment total or partial edentulous.

Bio Implant lines are used for both submerged technique is at the crest level, for Simple installation is strongly recommended submerged technique.

The Bio Implant dental implants should be considered an integral part of a complete implant system that includes in addition to Bio Implant accessories and the instruments also associated prosthetic components, for the V2 line characterized by unified platform are three different emergence profiles available: Narrow, Regular, Wide, not compatible Bio Implant dental implant Core V2 Ø2.9 which is expected instead dedicated prosthetic components.

Applications

The facilities must be selected according to the patient's condition and the operating technique. The latter is reflected in the surgical protocol, on the catalog. And a high degree of precision in the preparation of the implant site necessary through the respect of the sequence of the various steps provided by the surgical protocol. The thermal trauma from friction of the drills must be minimized with the use of sharp cutters, consequential, the appropriate number of revolutions, irrigation with sterile solutions to cool and healing times.

Contraindications

Insufficient amount of bone tissue or gum or insufficient or inadequate bone quality. Allergy to Titanium or its alloys with respect to patient selection to observe the same precautions and contraindications for dental and general surgery. Some of them include: reduced blood clotting (anticoagulants to antiplatelet therapies or, hereditary or acquired coagulation disorders) infection and inflammation of the mouth and the surrounding areas such as periodontitis and gingivitis (acute and chronic infections of the affected area, or infection inflammation of the soft tissues, bone infections, osteomyelitis). More over severe metabolic diseases (decompensated diabetes or inadequately treated, the calcium metabolism disorders, treatment with steroids and other drugs that interfere in bone metabolism, therapy with immune-suppressants, chemotherapy and radiotherapy, active endocrinological diseases, bone diseases, anatomical proximity or interference with nerve routes. Also insufficient gingival soft tissue for recasting, occlusion and / or dental articulation incorrect, limited interocclusal space. Furthermore psychic disorders, pain syndromes, poor oral hygiene, poor patient motivation to this type of treatment. Relative contraindications apply to bruxism, allergic patients, patients with additions to alcohol, drugs, nicotine.

Warnings

Dental implants brand Bio Implant can only be used by Dentists Doctors, dentists and maxillofacial surgeons, properly trained in the use of implant technique Bio Implant. Failure to observe the surgical protocol could determine the implant failure.

Dental implants brand Bio Implant should be considered an integral part of a complete implant system. Kristal recommend the use of accessories, tools and original prosthetic components.

Dental implants must be stored in its original packaging in a dry place, away from heat sources and at room temperature until the time of placement. Dental implants should not be used after the expiration date on the package.

Do not use the product if the package is opened or damaged. The lack of labeling and instructions may result in the use of the product beyond the expiration date and may be a necessary traceability. To modify the specification of the product packaging configuration could compromise sterility and functionality of the device. Dental implants are sold in sterile and single-use packaging, the reuse of single-use products creates a potential risk of infection in the patient or in those who use it.

The manufacturer is not liable in case of re-sterilization performed by third parties. Reprocessing may compromise the structural integrity of the product or cause it to malfunction. The removal of the plant from its packaging needs of the ideal about the aseptic standards. Contamination of the device may result in serious infection of the patient.

Kristal recommended to draw up appropriate clinical documentation. To ensure traceability of each plant is necessary to store information such as part number and lot.

Kristal provides a special label that shows the main features of the dental implant inserted to be placed on the patient's medical record.

To obtain the replacement of the goods when the quality does not match the standards of the manufacturer, if the wrong order or if visibly damaged packaging, it is necessary to warn by telephone or by e-mail within 10 days of receipt of goods; making the return of the product in its original filling a regular transport document of the commodity and the authorization request document made that comes from your target agent.

The devices removed from the oral cavity for a biological or mechanical failure must be treated for their disposal in biohazard waste. In general, please refer to the applicable regulations locally.

Precautions

In order to facilitate surgical planning it is advisable to carry out a thorough preliminary investigation to treatment screening, panoramic radiographs and possible T.A.C.

Side effects

Are as yet unknown side effects of Titanium C.P. Despite the high degree of success of the systems in titanium, they are not excluded the risks of failure or implant failure. Potential causes of these failures could be related to incorrect surgical techniques and failure to observe absolute contraindications to implantology.

Supply

Dental implants Plus, K-Plus, Core, K-Core and Core V2 Ø2.9: nr. 1 sterile plant with mounting device and pre-assembled cover screw.
Dental implants Core V2 and K-Core V2: nr. 1 sterile implant with cover screw.

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