

LANCER™ INSTRUMENTS PUT THE POWER AT YOUR FINGERTIPS



Quality instruments for every budget

QUALITY INSTRUMENTS FOR EVERY NEED AND BUDGET

Whether you're watching every expense, looking for high quality, or demanding the very best, Lancer offers instruments that will get the job done and look good doing it.



LANCER INSTRUMENTS

Lancer instruments are crafted from surgical stainless steel. Their screw joint design allows for easy tip alignment and adjustments.

These instruments are hand-buffed to a smooth, mirror luster for corrosion resistance and are designed to repeatedly withstand all accepted methods of sterilization without losing their aesthetic qualities.

Lancer instruments are available in a wide variety of the most popular designs and configurations.



STEALTH™ INSTRUMENTS

Constructed from the highest grade surgical stainless steel and the most modern, corrosion resistant materials available, Stealth instruments are designed to withstand all approved sterilization procedures.

Each Stealth instrument is hand-buffed and finished by true craftsmen to ensure the finest quality and performance. Stealth's uniquely effective designs dramatically increase stability without sacrificing range of motion, enhancing control while reducing required force.

Tips and cutting edges are manufactured from hardened steel or tungsten carbide for increased durability and better cutting performance over a longer period of time.

5 YEAR WARRANTY

Lancer warrants its instruments to be free of defects in material and workmanship for five (5) years. Misuse, abuse, or failure to properly maintain or care for the instrument will void warranty. Warranty does not cover routine maintenance, sharpening, or reconditioning. Cutting instruments should be reconditioned every 6-9 months. In the event of a warranty claim, it is the sole option of Lancer to repair or replace the instrument.

WIRE CUTTERS

Lancer's wire cutters have a variety of maximum cutting capacities up to .022x.028. Certain wire cutters are available with a 15 degree offset for easier access.



All Stealth™ wire cutters include Tungsten Carbide tip inserts.

**Hard Wire
(Round)**
.020

*Suggested maximum
wire thickness for
all Wire Cutters.*

**Stealth™
Hard Wire Cutter**

652-006



**Stealth™ Premium
Hard Wire Cutter**

652-006HD



**Stealth™
Power Cutter**

652-109



Hard/Soft Wire
1.5mm

*Suggested maximum
wire thickness for
Stealth™ Power Cutter*



**Lancer
Hard Wire Cutter**

626-121



**Lancer
Hard Wire - 15 Degree**

626-1000HWB

DISTAL END SAFETY CUTTERS

Lancer offers distal end cutters with the options of flush cut and safety hold.

Stealth™ distal end cutters include Tungsten Carbide tip inserts.




**Stealth™
Distal End**

652-101



**Stealth™
Distal End - Slim**

652-001



**Stealth™
Distal End - Slim - Long**

652-002



**Stealth™ Premium
Distal End - Slim**

652-001HD



**Stealth™ Premium
Distal End - Slim - Long**

652-002HD





All Stealth™ distal end safety cutters include Tungsten Carbide tip inserts.

Stealth™ Distal End Safety Cutter
Flush Cut - Long
652-188



Distal end cutters are available with standard or long handles.

**Stealth™ - Distal End
Flush Cut - Long**

652-188



**Stealth™
Distal End - Flush Cut**

652-102




Rectangular
.022 x .027


Round
.020


Braided
.016

*Suggested maximum
wire thicknesses for all
Distal End Safety Cutters.*



**Lancer
Distal End**

626-122



**Lancer
Distal End - Long**

626-122L

PIN AND LIGATURE CUTTERS


Soft Wire
.015

*Suggested maximum
wire thickness for all
Pin & Ligature Cutters.*

**Stealth™
Pin & Lig Cutter**

652-003



**Stealth™
Pin & Lig Cutter - Mini**

652-004



Lancer's pin and ligature cutters are available with standard or mini sized tips, or with a 45 degree angle for easier access to the posterior of the mouth.

Stealth™ pin and ligature cutters include Tungsten Carbide tip inserts.

**Stealth™ Premium
Pin & Lig Cutter**

652-003HD



**Stealth™ Premium
Pin & Lig Cutter - Mini**

652-004HD



45°

**Stealth™ - Pin & Lig
Cutter - 45 Degree**

652-005





*Stealth™ Premium
Pin and Ligature Cutter - Mini — 652-004HD*


**Stealth™ - Pin & Lig
Cutter w/ Spring**

652-104



**Stealth™ - Pin & Lig
Cutter - Mini w/ Spring**

652-105



45°

**Lancer - Pin & Lig
Cutter - 45 Degree**

626-117

15°

**Lancer - Pin & Lig Cutter
15 Degree**

626-118

7°

**Lancer - Pin & Lig Cutter
7 Degree**

626-119



**Lancer - Pin & Lig Cutter
Micro**

626-120

GOOSE NECK PLIERS



**Stealth™
Goose Neck Pliers**

652-184

WEINGART PLIERS

Our versatile Weingart pliers are primarily used for insertion and removal of archwires. Options of straight, angled or tapered tips are available, as well as tungsten carbide.




**Lancer
Weingart - Angled**

626-124



**Stealth™
Weingart - Regular**

652-117



**Stealth™
Weingart Pliers w/
Tungsten Carbide Tips**

652-118

TC Tungsten
Carbide

HOW PLIERS

*Stealth™ How Plier
Regular — 652-113*



**Stealth™ - How Plier
Titanium Alloy**

652-112



**Stealth™
How Plier - Regular**

652-113

Our How pliers feature serrated tips to facilitate the placement and removal of archwires, pins and other auxiliary products.

Offered with options of straight tips, offset tips and in Titanium alloy.

**Stealth™
How Plier - Offset Tips**

652-115




**Lancer
How Plier - Small Tips**

626-123



**Lancer
How Plier - Offset Tips**


626-140



BIRD BEAK PLIERS

These versatile loop forming pliers allow the clinician to bend wires without the risk of damage.

Bird Beak pliers are available in a variety of tip length, diameters and wire bending capacities.




**Stealth™
Bird Beak**

652-009



**Stealth™
Bird Beak w/ Cutter**


652-120



**Stealth™ - Bird Beak
w/ Tungsten Carbide Tips**

652-121

MTC Tungsten
Carbide



**Lancer
Bird Beak**

626-129



**Lancer - Bird Beak w/
Square Serrated Tips**

631-023

**Max Bending
Capacity**
.030

*Suggested maximum
wire thickness for
all Bird Beak Pliers.*

FORMING PLIERS

Hollow Chop pliers are used for subtle contouring of the archwire up to size .030 and forming palatal arches.

Lingual arch forming pliers are developed to form double and triple bends in the lingual archwire.



**Stealth™
Hollow Chop Plier**

652-014



**Lancer
Hollow Chop Plier**

626-126



**Stealth™ Arch Forming
Plier - Grooved**

652-149



**Lancer
Lingual Arch Forming Plier**

626-130

LIGHT WIRE PLIERS

These versatile loop forming pliers allow the clinician to bend wires without the risk of damage.

**Stealth™ - Light Wire
w/ Grooved Square Tip**

652-124



**Stealth™ - Light Wire
w/ Tungsten Carbide Tips**

652-126

WIRE BENDING PLIERS



**Stealth™
Nance Loop Closing Plier**

652-136

Nance pliers allow the clinician to form precise loops of various sizes in rectangular wire.



**Lancer
Tweed Arch Bending Plier**

626-136



**Stealth™
Tweed Arch Bending Plier**

652-012



**Stealth™
Adams Plier**

652-141




**Stealth™
Wire Bending Plier**

652-144



**Stealth™
Stop Plier**

652-142



Lancer - Stop Plier

626-103



**Stealth™
Heavy Three Jaw Plier**

652-013



**Stealth™
Three Jaw Plier**

652-139



**Lancer
Three Jaw Plier**

626-125

Lancer's three-jaw pliers are designed for contouring and forming all archwires, especially NiTi wire.

LOOP FORMING PLIERS

These pliers are ideal for forming helical loops and intraoral tip back bends. Some of these pliers have removable replacement tips. The maximum bending capacity ranges from .020 to .022x.028.



**Stealth™
Tweed Loop Forming**

652-129



**Stealth™ - Loop Forming
w/ Replaceable Tip**

652-011



**Stealth™
Loop Forming**

652-131



**Stealth™
Jaraback Plier**

652-128

**Stealth™
2-Step Replacement Tip
(2 tips, 2 screws, hex key)**

652-135

**Lancer
Tweed Loop Forming**

626-127

**Lancer - Tweed
3-Step Replacement Tip
(2 tips, 2 screws, hex key)**

626-352T



**Lancer
Optical Plier**

631-301

LIGATING INSTRUMENTS

The Mathieu and Mosquito hemostats feature tapered serrated tips and locking handles especially designed for the ligation process. These features allow the clinician to firmly grip ligatures without intruding into the lumen.



**Stealth™
Ligature Director**


652-125



**Stealth™
Ligature Tying Plier
Coon Style**

652-023

Coon Style ligature tying pliers are designed to automatically lock wire ends as tips spread when tying preformed steel ligatures.




**Lancer
Mosquito Hemostat**

626-306

**Lancer Mathieu - Ultra
Fine Tip - Straight**

626-138



**Stealth™
Mathieu Hemostat**

652-024

**Stealth™ - Elastic
Placing Plier - Hook Tip**


652-154

**Stealth™
Ligating Plier**

652-155

DEBONDING INSTRUMENTS

Lancer's bracket removal pliers are specifically designed to debond brackets with minimal force and patient discomfort. These pliers are available with straight or angled tips for optimal ergonomics.



**Lancer - Bracket
Removing Plier**

626-135

**Lancer
Bracket Removing Plier
Replacement Blade 2/pk**

626-135T

**Lancer - Bracket
Removing Plier - Straight**

626-346



**Stealth™ - Direct Bond
Removing Plier - Angled**

652-019



**Stealth™ - Direct Bond
Removing Plier - Straight**

652-020



**Stealth™
Bond Removing Plier**

652-021



**Stealth™
Adhesive Removing Plier**

652-183



**Stealth™ - Adhesive
Removing Replacement Tip**

652-185

BAND TOOLS

The tips and occlusal pads of Lancer's band removal pliers are accurately aligned to facilitate easy removal of bands in the posterior of the mouth. Available with long or short tips and with replaceable pads.

Lancer - Band Seater / Instrument

626-312

Lancer - Band Removing Plier - Short

626-133

Band Removing Plier Replacement Pads 4.66mm 3/pk

626-134T

Lancer - Band Removing Plier - Long

626-134

Stealth™ - Posterior Band Removing Plier

652-025

Lancer - Bite Stick Heat Sterilizable 1/pk

631-111

Lancer Band Pusher

Blue 631-114

Green 631-115

Lancer Band Pusher

631-307

Stealth™ Premium Posterior Band Removing Plier

652-025HD

Stealth™ - Posterior Band Removing Plier Long

652-172

BRACKET TOOLS



**Stealth™
Boone Bracket
Positioning Gauge**

652-152



**Lancer - Deluxe
Bracket Holder**

626-102



**Lancer
Bracketizer**


626-104

**Lancer - Bracket
Positioning Gauge .018**

631-316

**Lancer - Bracket
Positioning Gauge .022**

631-317



**Lancer
Gram Gauge
Dontrix 500gm**

631-320



TWEEZERS

**Lancer - Bracket
Placing - Chrome**

626-100



**Lancer
Keat Tweezer**

626-116



**Lancer - Bracket
Placing - Satin**

626-309



**Lancer - Bracket
Placing - Slim**

631-330



**Stealth™ - Bracket
Placing Tweezer**

652-018



CRIMPING PLIERS



**Stealth™
Hook Crimping Plier**


652-017



**Lancer
Hook Crimping Plier**


626-115

OTHER INSTRUMENTS



**Lancer
Torque Plier**

626-131



**Lancer
Finishing Plier 3/4 mm**

626-150




**Stealth™ - Self-Ligating
Bracket Instrument**

652-028

**Stealth™
Teflon Pad
Replacement 3/pk**

652-175



**Stealth™
Separating Plier**

652-182

Tie back pliers provide the clinician with an efficient technique to permanently bend NiTi wire distal to the molar tubes.




**Lancer
1mm Step Plier**

626-101



**Lancer - Explorer /
Elastic Remover**

463-998



**Stealth™
Tie Back Plier**

652-186



ORTHOBOX™

Lancer's OrthoBox™ provides a safe and efficient way to handle, carry and sterilize instruments.

- Manufactured of high-quality stainless steel with a rust-resistant anodized aluminum cover
- Can be inserted directly into autoclave
- Holds up to 14 instruments
- L 9.85 in x W 5.5 in H 2.4 in, (L 25 cm x W 14 cm x H 6 cm)

**OrthoBox™
Blue**

626-323

**OrthoBox™
Red**

626-324

**OrthoBox™
Green**

626-325

**OrthoBox™
Yellow**

626-326

Colors are for demonstration only. Colors are applied to anodized aluminum and will not appear exactly as shown.



(800) 854-2896 • www.lancerortho.com

International: +1 (760) 744-5585 • Fax: +1 (760) 598-0418
2726 Loker Avenue West • Carlsbad, California 92010 (U.S.A.)
info@lancerortho.com