

# SPLINX™ RETAINER WIRE MAKES PERMANENT RETENTION EASY.



- Suitable for all types of retention and splinting both upper and lower arches
- Made from a flattened woven wire coaxial (6 wires) for an easier adherence to tooth surface without slipping
- Completely passive, easily shaped and with no memory residual strength
- Compatible with all adhesives
- Wire Length: 8" (20.4 cm); Wire Thickness: 0.25 mm; Wire Height: 0.69 mm
- Tube of 8 wires (1.63 meters total)

*The first bonded lingual wire that doesn't fray when cut.*



SPLINX™  
RETAINER WIRE

Tube of 8  
8" Lengths

0.0155

550-120



(800) 854-2896 • [www.lancerortho.com](http://www.lancerortho.com)

International: +1 (760) 744-5585 · Fax: +1 (760) 598-0418  
2726 Loker Avenue West · Carlsbad, California 92010 (U.S.A.)  
[info@lancerortho.com](mailto:info@lancerortho.com)