

VELOCITY™ BUCCAL TUBES METICULOUSLY ENGINEERED FOR PERFORMANCE



*Some of the most trusted
buccal tubes on the market.*

Velocity™ Buccal Tubes

Manufactured using the latest Metal Injection Molding (MIM) technologies, Velocity™ buccal tubes offer a stylish, low-profile design for maximum patient comfort and clinical effectiveness.

Lancer's Velocity™ line of bondable and weldable convertible and non-convertible buccal tubes incorporates the latest technology in design and manufacturing. Our convertible tubes feature an integrated molded cap to ensure conversion only when clinically removed.

Velocity buccal tubes answer the growing demand for comfortable, easy-to-use appliances, resulting in efficient treatment and reduced chair time.



VELOCITY™ SINGLE BUCCAL TUBES

The latest design of an extremely low-profile buccal tube, Velocity™ buccal tubes are setting the standard in exceeding clinical demands and assuring patient comfort. The funneled entrance of our Velocity™ tubes facilitates easy wire insertion. Wires are guided effortlessly into the broad entrance, simplifying arch wire adjustments and saving valuable chair time.

This tube offers patient comfort and clinical effectiveness in one powerful appliance.

- ***Funneled entrance for easy wire insertion, occlusal/gingival positioning guides to aid in placement, and a distal notch for additional ligation***
- ***Bondable and Weldable – Available in Roth, MBT and Standard Rx's, .018 and .022***
- ***Tube length 4.3mm***

VELOCITY™ SINGLE BUCCAL TUBES

1ST / 2ND MOLAR

Upper / Maxillary

				Long Bondable		Short Bondable	
				.018 Each	.022 Each	.018 Each	.022 Each
Torque	Offset	L/R		3DL8-SN05	3DL2-SN05	3DS8-SN05	3DS2-SN05
				3DL8-SN06	3DL2-SN06	3DS8-SN06	3DS2-SN06
Roth	-10°	8°	R	3DL8-SN01	3DL2-SN01	3DS8-SN01	3DS2-SN01
			L	3DL8-SN02	3DL2-SN02	3DS8-SN02	3DS2-SN02
MBT	-14°	10°	R				
			L				

Lower / Mandibular

				Long Bondable		Short Bondable	
				.018 Each	.022 Each	.018 Each	.022 Each
Torque	Offset	L/R		3DL8-SN07	3DL2-SN07	3DS8-SN07	3DS2-SN07
				3DL8-SN08	3DL2-SN08	3DS8-SN08	3DS2-SN08
Roth	-25°	5°	R	3DL8-SN03	3DL2-SN03	3DS8-SN03	3DS2-SN03
			L	3DL8-SN04	3DL2-SN04	3DS8-SN04	3DS2-SN04
MBT	-20°	0°	R				
			L				

EDGEWISE/STANDARD 1ST / 2ND MOLAR

Upper / Lower

				Long Bondable		Short Bondable	
				.018 Each	.022 Each	.018 Each	.022 Each
Torque	Offset	L/R		3DL8-SN09	3DL2-SN09	3DS8-SN09	3DS2-SN09
				3DL8-SN10	3DL2-SN10	3DS8-SN10	3DS2-SN10
Edgewise/ Standard	0°	0°	UR/LL				
			UL/LR				

WELDABLE BUCCAL TUBES



Upper / Maxillary

				.018 Each	.022 Each
				Torque	Offset
Std	0°	0°	R	3W8-SN09	3W2-SN09
			L	3W8-SN10	3W2-SN10
Roth	-10°	8°	R	3W8-SN05	3W2-SN05
			L	3W8-SN06	3W2-SN06
MBT	-14°	8°	R	3W8-SN01	3W2-SN01
			L	3W8-SN02	3W2-SN02



Lower / Mandibular

				.018 Each	.022 Each
				Torque	Offset
Std	0°	0°	R	3W8-SN10	3W2-SN10
			L	3W8-SN09	3W2-SN09
Roth	-25°	5°	R	3W8-SN07	3W2-SN07
			L	3W8-SN08	3W2-SN08
MBT	-20°	0°	R	3W8-SN03	3W2-SN03
			L	3W8-SN04	3W2-SN04



VELOCITY™ SINGLE CONVERTIBLE BUCCAL TUBES

Produced from state-of-the-art design and MIM manufacturing technologies, Lancer's Velocity Single Convertible Buccal Tubes feature an integrated molded cap to ensure conversion only when clinically removed. Following archwire slot conversion, the twin tie wing design facilitates secure ligation.

- **Bondable and Weldable – Available in Roth and MBT Rx's, .018 and .022 slot sizes**
- **Reduced occlusal interference and gingival impingement**
- **Mesial extension tab facilitates removal of the "peel-away" cap**
- **Tube length: 3.5mm**

VELOCITY™ SINGLE CONVERTIBLE BUCCAL TUBES

Upper / Maxillary

	Torque	Angle	L/R	Long Bondable		Short Bondable	
				.018 Each	.022 Each	.018 Each	.022 Each
Roth	-10°	8°	R	3DL8-SC05	3DL2-SC05	3DS8-SC05	3DS2-SC05
			L	3DL8-SC06	3DL2-SC06	3DS8-SC06	3DS2-SC06
MBT	-14°	10°	R	3DL8-SC01	3DL2-SC01	3DS8-SC01	3DS2-SC01
			L	3DL8-SC02	3DL2-SC02	3DS8-SC02	3DS2-SC02

Lower / Mandibular

Roth	-25°	5°	R	3DL8-SC07	3DL2-SC07	3DS8-SC07	3DS2-SC07
			L	3DL8-SC08	3DL2-SC08	3DS8-SC08	3DS2-SC08
MBT	-20°	0	R	3DL8-SC03	3DL2-SC03	3DS8-SC03	3DS2-SC03
			L	3DL8-SC04	3DL2-SC04	3DS8-SC04	3DS2-SC04

WELDABLE BUCCAL TUBES

Upper / Maxillary

	Torque	Offset	L/R	.018 Each	.022 Each
Roth	-10°	8°	R	3W8-SC05	3W2-SC05
			L	3W8-SC06	3W2-SC06
MBT	-14°	10°	R	3W8-SC01	3W2-SC01
			L	3W8-SC02	3W2-SC02

Lower / Mandibular

	Torque	Offset	L/R	.018 Each	.022 Each
Roth	-25°	5°	R	3W8-SC07	3W2-SC07
			L	3W8-SC08	3W2-SC08
MBT	-20°	0	R	3W8-SC03	3W2-SC03
			L	3W8-SC04	3W2-SC04

VELOCITY™ DOUBLE CONVERTIBLE BUCCAL TUBES

Latest design. MIM technology. Integrated molded cap ensures conversion only when clinically removed.

- **Bondable and Weldable – Available in Roth and MBT Rx's, .018 and .022 with auxiliary .018 x .025 tube**

Incorporating a twin tie wing design for secure ligation following archwire slot conversion, Lancer's Velocity Double Convertible Buccal tubes are also offset from the gingival to avoid tissue impingement.

Velocity™ Double Convertible Buccal tubes feature a .018 x .025 auxiliary arch wire slot used for piggyback techniques.



VELOCITY™ DOUBLE CONVERTIBLE BUCCAL TUBES

1ST / 2ND MOLAR

Upper / Maxillary

	Torque	Offset	Aux	L/R	.018 Each	.022 Each
Roth	-10°	8°	.018 x .025	R	3DL8-DC05	3DL2-DC05
				L	3DL8-DC06	3DL2-DC06
MBT	-14°	10°	.018 x .025	R	3DL8-DC01	3DL2-DC01
				L	3DL8-DC02	3DL2-DC02

Lower / Mandibular

	Torque	Offset	Aux	L/R	.018 Each	.022 Each
Roth	-25°	5°	.018 x .025	R	3DL8-DC07	3DL2-DC07
				L	3DL8-DC08	3DL2-DC08

WELDABLE BUCCAL TUBES

Upper / Maxillary

	Torque	Offset	Aux	L/R	.018 Each	.022 Each
Roth	-10°	8°	.018 x .025	R	3W8-DC05	3W2-DC05
				L	3W8-DC06	3W2-DC06
MBT	-14°	10°	.018 x .025	R	3W8-DC01	3W2-DC01
				L	3W8-DC02	3W2-DC02

Lower / Mandibular

	Torque	Offset	L/R	.018 Each	.022 Each
Roth	-25°	5°	R	3W8-DC07	3W2-DC07
			L	3W8-DC08	3W2-DC08



VELOCITY™ MINI BUCCAL TUBES

Velocity Mini Buccal Tubes are small and powerful, incorporating all the functional elements of our standard Velocity series into an attractive, compact design that allows for earlier bonding of partially erupted second molars.

Designed with a funneled entrance to facilitate easy wire insertion, Velocity Mini Buccal Tubes allow wires to be guided gently and easily into the entrance.

Bonds easily by the clinician with less trauma to the soft tissue.

- Fits easily and precisely on the mesial buccal cusp of 2nd molars
- Funneled entrance for easy wire insertion, occlusal/gingival positioning guides to aid in placement
- Smooth, polished comfortable hook for ultimate patient well-being
- 80 gauge mesh, compound contour pad provide more surface area for greater bond strength
- Tube length 3mm

VELOCITY™ MINI BUCCAL TUBES

Upper / Maxillary

	Torque	Offset	L/R	.018 Each	.022 Each
Roth	-10°	0°	R	3DS8-SN15	3DS2-SN15
			L	3DS8-SN16	3DS2-SN16
MBT	-14°	0°	R	3DS8-SN11	3DS2-SN11
			L	3DS8-SN12	3DS2-SN12

Lower / Mandibular

	Torque	Offset	L/R	.018 Each	.022 Each
Roth	-25°	0°	R	3DS8-SN17	3DS2-SN17
			L	3DS8-SN18	3DS2-SN18
MBT	-20°	0°	R	3DS8-SN13	3DS2-SN13
			L	3DS8-SN14	3DS2-SN14

Upper/Lower

	Torque	Offset	L/R	.018 Each	.022 Each
Edgewise/ Standard	0°	0°	UR/LL	3DS8-SN19	3DS2-SN19
			UL/LR	3DS8-SN20	3DS2-SN20



VELOCITY™ TRIPLE CONVERTIBLE BUCCAL TUBES

Latest design. MIM technology. Integrated molded cap ensures conversion only when clinically removed.

Velocity triple tubes feature a single tube along with an auxiliary 18x25 arch wire slot and a .045 or .051 headgear tube. Our occlusal tubes are ideal for cervical headgear or appliances that need center of resistance on the occlusal plane.

- **3.5mm M/D – Weldable – Available in Roth and MBT Rx's, .018 and .022**

VELOCITY™ TRIPLE CONVERTIBLE BUCCAL TUBES

These superior-performing buccal tubes feature a low-profile design for reduced occlusal interference, comfortable hooks for maximum patient comfort, and a practical base design that anatomically conforms to the tooth. Their convenient “peel-off” convertible caps are integrated to avoid unintentional conversion during treatment.

WELDABLE BUCCAL TUBES

1ST MOLAR

Upper / Maxillary

					.018	.022
					Each	Each
Roth	Torque	Offset	HGO	L/R	3W8-TC05	3W2-TC05
				L	3W8-TC06	3W2-TC06
	-10°	8°	0.045°	R	3W8-TC15	3W2-TC15
				L	3W8-TC16	3W2-TC16
MBT	-14°	10°	0.045°	R	3W8-TC01	3W2-TC01
				L	3W8-TC02	3W2-TC02
	-14°	10°	0.051°	R	3W8-TC11	3W2-TC11
				L	3W8-TC12	3W2-TC12

LANCER PRODUCES SOME OF THE MOST POPULAR AND TRUSTED BUCCAL TUBES ON THE MARKET.

With buccal tubes being such a critical component in orthodontic treatment, Lancer has meticulously engineered a full complement of weldable and bondable convertible and non-convertible tubes for all your treatment needs. These buccal tubes are available both individually, and custom-welded directly onto bands to meet your specific clinical needs.

As with all the appliances we manufacture, Lancer applies internationally recognized quality standards in the production of our tubes to help ensure that we consistently deliver world-class quality to our customers.



(800) 854-2896 • www.lancerortho.com

International: +1 (760) 744-5585 • Fax: +1 (760) 598-0418
2726 Loker Avenue West • Carlsbad, California 92010 (U.S.A.)
info@lancerortho.com